

AIRWAY ASSISTANCE

NASOPHARYNGEAL AIRWAY

Nasopharyngeal airways provides flexible and simple access to secure a patient's airway. As the name implies, it is inserted through the nose and can be used on both conscious and unconscious patients. The tube has a flared proximal end to prevent it from entering too deeply. The 170 mm length allows it to reach the back of the throat providing clear air passage. Each tube has an adjustable flange to adjust the length. To ease in insertion each Nasopharyngeal tubes come with a 5g package of lubricant. Three models are provided: Silkolatex, Latex-free and soft rubber. The more common sizes for adult use are the 26 and 28 GR, but other sizes are available.



NASOPHARYNGEAL SILKOLATEX

Adjustable Length
Single-Use
Length 170 mm



NASOPHARYNGEAL LATEX-FREE

Adjustable Length
Single-Use
Length 170 mm



NASOPHARYNGEAL SOFT RUBBER

Adjustable Length
Re-Usable
Length 170 mm

Article Nr.	Military NSN	Description	Unit	Price (excl. VAT)
185200 GR. 26	-	Nasopharyngeal Silkolatex 26F - 170 mm	1	€ 6.74
185200 GR. 28	-	Nasopharyngeal Silkolatex 28F - 170 mm	1	€ 6.74
185420 GR. 26	-	Nasopharyngeal Latex-Free 26F - 170 mm	1	€ 8.72
185420 GR. 28	-	Nasopharyngeal Latex-Free 28F - 170 mm	1	€ 8.72
125200 GR. 26	-	Nasopharyngeal Soft Rubber 26F - 170 mm	1	€ 14.25
125200 GR. 28	-	Nasopharyngeal Soft Rubber 28F - 170 mm	1	€ 14.25

MEDICAL SALES CONSULTANTS AUSTRIA GMBH

Wachaustrasse 13, 3123 Neustift, Austria
T: +43 2786 631 38, E-Mail: info@medicalsei.com
Web: www.medicalscia.at

MEDICAL SALES CONSULTANTS GERMANY GMBH

Vor dem Schonekindtor 13, 59494 Soest, Germany
T: +49 2921 969 78 50, E-Mail: info@medicalseg.de
Web: www.medicalscg.de