

TACTICAL RESCUE


SMART XTRACT SR COMPLETE & ELEMENT PROTECTION SYSTEM

The XTract® SR is a lightweight, compact rescue stretcher designed for rapid extraction of a seriously injured casualty from a hostile environment. It provides a modular approach, providing the rescuer the ability to configure and use as required for each situation. It is designed to cocoon the casualty, helping to secure and protect medical interventions and to prevent secondary injuries.

The XTract® SR "Complete" System comes in a carrying bag which includes the XTract® SR stretcher with a sled base for drag extractions and straps for lifting, dragging or hands-free carrying. An "Element Protection System" (EPS) is available as a separate item that can be used as a cover over the casualty for added warmth and water protection.



XTRACT® SR COMPLETE SYSTEM
 DIMENSIONS AND WEIGHT
 Packed: 56 x 26 cm, Weight: 4.5 kg
 Xtract SR Opened: 96 x 208 cm
 Max. Carrying Load: 300 kg



XTRACT® SR - Weight: 1.25 kg



EPS - Weight: 400 g



XTRACT® SR Complete System

Article Nr.	Military NSN	Description	Unit	Price (excl. VAT)
72111	-	XTract® SR Stretcher - Black	1	€ 364.80
72121	-	XTract® SR Complete System - Black	1	€ 747.20
72116	-	Element Protection System (EPS) - Green	1	€ 142.40