

TACTICAL RESCUE SAVIOUR TECHNICAL STRETCHER

The Saviour Technical Stretcher is a rapid response stretcher to evacuate casualties safely from hazardous areas. Up to six persons can carry the Saviour and there are eyelets that can be used with the lifting harness or drag straps. This allows for vertical lifting and lowering, as well as dragging or hauling along evacuation routes. The Saviour Technical may also be lifted horizontally using an approved lifting strap. The Saviour stretcher is made of low-density polyethylene, is X-Ray translucent and can also be used in CT scanners.

FEATURES AND ADVANTAGES

- > Works well in wet or sandy environments
- > Floats when dropped overboard, unloaded
- > Floats when loaded with a patient in a lifejacket
- > Features chest straps that cross the thorax and a figure 8 harness for vertical lifting
- > 6 soft carrying handles
- > Can be carried by 2 or 4 people and if needed carrying poles can be threaded through the handles when only 2 rescuers are available
- > Colored coded straps for fast and easy connection
- > All straps can be easily replaced
- > After patient use, rinse the stretcher with cold water and if needed a low pressure washer can be used
- > Body fluids should be washed from the stretcher and straps with a non-corrosive disinfectant and then rinsed in fresh water



SAVIOUR TECHNICAL STRETCHER

Dimensions Rolled: 70 x 30 cm
 Opened: 215 x 63 cm, Weight incl. Bag: 7 kg
 Weight w/o Bag: 5.7 kg
 Max. Carrying Load: 200 kg



SPINAL SPLINT

Opened: 102 x 20 x 44cm Weight: 3.5kg
 Max. Carrying Load: 200kg

Article Nr.	Military NSN	Description	Unit	Price (excl. VAT)
65263	-	Saviour Technical Stretcher	1	€ 750.00
66532	-	Saviour Technical Stretcher - Adjustable Strap for Hoisting	1	€ 226.00
SVDS	-	Saviour Technical Stretcher - Drag Strap	1	€ 30.00
83980-0	-	Spinal Splint for Saviour Technical Stretcher	1	€ 480.00
83985-0	-	Spinal Splint Carrying Bag	1	€ 47.70

MEDICAL SALES CONSULTANTS AUSTRIA GMBH

Wachaustrasse 13, 3123 Neustift, Austria
 T: +43 2786 631 38, E-Mail: info@medicalseci.com
 Web: www.medicalscia.at

MEDICAL SALES CONSULTANTS GERMANY GMBH

Vor dem Schonekindtor 13, 59494 Soest, Germany
 T: +49 2921 969 78 50, E-Mail: info@medicalsecg.de
 Web: www.medicalscg.de