

HEMOSTASIS AGENTS

QUIKLOT® EMS ROLLED GAUZE™ QUIKLOT® EMS DRESSING

QuikClot® EMS products are the same as all other well-known QuikClot® kaolin products, but the difference is they are specifically packaged for civilian pre-hospital use. There is no X-Ray detectable strip and the package has a tear-open top allowing for the fastest access to the product, which can be especially important if the user is wearing thick gloves. As with all QuikClot® products, the kaolin activates the patients coagulation factors to accelerate the time to "stop the bleeding".

Unlike other hemostatic agents, QuikClot® products contain no animal components, alleviating concerns of possible allergic reactions. And like all QuikClot® products it conforms to the wound. QuikClot® EMS products can be used superficially or for deep wounds. All QuikClot® products are effective on patients who take anticoagulant medication.



QUIKLOT® EMS DRESSING

Opened: 10 x 10 cm

QUIKLOT® EMS ROLLED GAUZE™

Opened: 7.6 x 122 cm

Article Nr.	Military NSN	Description	Unit	Price (excl. VAT)
491	-	QuikClot® EMS Dressing	Box 3	€ 24.99
493	-	QuikClot® EMS Rolled Gauze™	1	€ 16.67